Technovit® EPOX

Technovit® EPOX

- D** Gebrauchsinformation
- GB** User Instructions
- F** Mode d'emploi
- NL** Gebruiksaanwijzing
- E** Modo de empleo
- I** Istruzioni per l'uso
- RU** Инструкция по применению
- CN** 用户说明书

Heraeus

Heraeus Kulzer GmbH
 Philipp-Reis-Straße 8/13
 D-61273 Wehrheim/Ts.
 Telefon (0 60 81) 9 59-0
 Telefax (0 60 81) 9 59-398
 e-mail: technik.wehrheim@heraeus.com
 Internet: www.Kulzer-Technik.com

Kulzer:

Stand: März 2009
66040445/00 DOG/9-spr/03.09

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Art.-Nr. 66040438 –
 Technovit EPOX Hardener regular,
 1 x 500 g

Art.-Nr. 66040439 –
 Technovit EPOX Hardener fast,
 1 x 500 g

GB User Instructions Technovit EPOX

Technovit EPOX is a Epoxy-Resin-System, consisting of a resins component „**Technovit EPOX Resin**“ and two different hardener components „**Technovit EPOX Hardener fast**“ or „**Technovit EPOX Hardener regular**“.

Area of application is embedding of materialographic samples, especially porous materials, which has to be embedded under vacuum. Technovit EPOX characterizes itself through high transparency, strong adhesion on the sample as well as low shrinkage and UV-stability.

Sample preparation

Clean samples carefully! The sample surface has to be free of dust, grease or other contamination. Otherwise curing can be disturbed or bubbles appear on the sample surface.

Measuring and mixing of the components

In order to achieve an optimal result, the mixing ratio has to be adhered exactly. The ratio is two parts of „Resin“ to one part of „Hardener“ by weight.

Place „Technovit EPOX Resin“ in a mixing jar. Afterwards weigh in the „Technovit EPOX Hardener“. The composition has to be mixed very until the material is free of streaks (stir minimum for 60 sec.). Don't stir too heavy in order to avoid inclusion of bubbles in the material.

After mixing, let set down the material for round about 60 sec. to allow leaking the bubbles. Now pour slowly the material over the specimen, prepared in the embedding mould. Take care that there remains no air under the bottom of the specimen.

Quantities and curing times refer to separate chart!

Storage

Store the materials in a well ventilated place. Do not store at temperatures exceeding 25 °C. Keep away from children. Do not use after the expiry date.

Danger-/Safetey advice

Technovit EPOX Resin contains:
 Bisphenole A Epoxyresin, liquid (MG < 700),
 Bisphenole-F Epichlorhydrinresin (MG < 700),
 Oxirane, Mono
 [(C12-14-alkyloxy)methyl]derivate
 Irritating to respiratory system and skin. May cause sensitization by skin contact. Toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment. Do not breathe vapour. Avoid contact with the skin. In case of contact with eyes, rinse immediately with plenty of water and seek medical

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advice. After contact with skin, wash immediately with plenty of water. Wear suitable gloves and eye/face protection. Use appropriate container to avoid environmental contamination.

Technovit EPOX Hardener contains:
 3-Aminomethyl-3,5,5-trimethylcyclohexylamine, Tricyclododecandiamine, Trimethylhexan-1,6-diamin
 Benzylalkohole

Harmful by inhalation, in contact with skin and if swallowed. Causes burns. May cause sensitization by skin contact. Harmful to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

Keep container in a well-ventilated place. When using do not eat or drink. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. Wear suitable protective clothing, gloves and eye/face protection. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible). Avoid release to the environment. Refer to special instructions/safety data sheets.

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Art. no. 66040437 –
 Technovit EPOX Resin,
 1 x 1000 g

Art. no. 66040438 –
 Technovit EPOX Hardener regular,
 1 x 500 g

Art. no. 66040439 –
 Technovit EPOX Hardener fast,
 1 x 500 g

Aushärungsdaten/Curing datas								
Einbettform/ Embedding mould	Menge ohne Probe/Quantity without Specimen	Umgebungs- Temperatur/ Environment temperature	Technovit EPOX Resin/ Technovit EPOX hardener regular			Technovit EPOX Resin/ Technovit EPOX hardener fast		
			Max. Aushärte- Temperatur/ max. curing temperature (Tmax)	Zeit bis/ Time until Tmax	Zeit bis Ende messbare Temperatur- Entwicklung/ Time until measurable temperature development	Max. Aushärte- Temperatur/ max. curing temperature (Tmax)	Zeit bis/ Time until Tmax	Zeit bis Ende messbare Temperatur- Entwicklung/ Time until measurable temperature development
25 mm	12 g (8 g Resin/4 g Hardener)	20 °C	30 °C	120 Min.	ca. 18 h	37 °C	90 Min.	ca. 10 h
30 mm	18 g (12 g Resin/6 g Hardener)	20 °C	35 °C	110 Min.	ca. 18 h	57 °C	80 Min.	ca. 10 h
40 mm	30 g (20 g Resin/10 g Hardener)	20 °C	45 °C	105 Min.	ca. 18 h	110 °C	70 Min.	ca. 10 h
50 mm	45 g (30 g Resin/15 g Hardener)	20 °C	65 °C	100 Min.	ca. 18 h	144 °C	60 Min.	ca. 10 h